## PED, NEO, PDL, PDC, PLC

English	4
Français	8
Deutsch	13
Italiano	18
Español	22
Português	27
Nederlands	31
Svenska	35
Dansk	39
Norsk	43
Suomi	47
Ελληνικά	51
Symbols	56

Table 1.

	PED	NEO	PDL	PDC	PLC
	X	X	X	X	X
	X	X	X	X	X

Table 2. PED & PDL



			
Product Size	I.D. (mm)	O.D. (mm)	Length (mm)
3.0 PED	3.0	4.5	39
3.5 PED	3.5	5.2	40
4.0 PED	4.0	5.9	41
4.5 PED	4.5	6.5	42
5.0 PED	5.0	7.1	44
5.5 PED	5.5	7.7	46
5.0 PDL	5.0	7.1	50
5.5 PDL	5.5	7.7	52
6.0 PDL	6.0	8.3	54
6.5 PDL	6.5	9.0	56


Table 3. NEO

			
Product Size	I.D. (mm)	O.D. (mm)	Length (mm)
3.0 NEO	3.0	4.5	30
3.5 NEO	3.5	5.2	32
4.0 NEO	4.0	5.9	34
4.5 NEO	4.5	6.5	36

**Table 4. PDC & PLC**

					
Product Size	I.D. (mm)	O.D. (mm)	Length (mm)	Leak Test Inflation Volume (cc)	Cuff Resting Diameter (mm)
4.0 PDC	4.0	5.9	41	8.0	12.4
4.5 PDC	4.5	6.5	42	8.0	12.5
5.0 PDC	5.0	7.1	44	9.0	15.0
5.5 PDC	5.5	7.7	46	9.0	14.9
5.0 PLC	5.0	7.1	50	9.0	15.0
5.5 PLC	5.5	7.7	52	9.0	14.9
6.0 PLC	6.0	8.3	54	11.0	19.3
6.5 PLC	6.5	9.0	56	11.0	19.3

**Table 5.**

			
Product Size	A (mm)	Product Size	A (mm)
3.0 NEO	9.6	4.5 PED & PDC	13.9
3.5 NEO	10.3	5.0 PED & PDC	14.8
4.0 NEO	11.1	5.5 PED & PDC	15.8
4.5 NEO	11.9	5.0 PDL & PLC	20.8
3.0 PED	11.4	5.5 PDL & PLC	21.8
3.5 PED	12.9	6.0 PDL & PLC	22.8
4.0 PED & PDC	13.4	6.5 PDL & PLC	23.7

This instruction insert applies to the following Shiley Tracheostomy products:

PED	Pediatric Tracheostomy Tube
NEO	Neonatal Tracheostomy Tube
PDL	Long Pediatric Tracheostomy Tube
PDC	Cuffed Pediatric Tracheostomy Tube
PLC	Cuffed Long Pediatric Tracheostomy Tube

#### Instructions for Use

##### **WARNINGS:**

CAREFULLY READ THESE INSTRUCTIONS PRIOR TO USE. THESE INSTRUCTIONS ARE APPLICABLE ONLY TO THE SHILEY PEDIATRIC TRACHEOSTOMY TUBE (PED), SHILEY LONG PEDIATRIC TRACHEOSTOMY TUBE (PDL), SHILEY NEONATAL TRACHEOSTOMY TUBE (NEO), SHILEY PEDIATRIC WITH CUFF (PDC), AND SHILEY LONG PEDIATRIC WITH CUFF (PLC) AND SHOULD BE READ CAREFULLY BEFORE USING THE TUBE.

STERILE ONLY IF THE PROTECTIVE PACKAGE IS NOT OPEN, DAMAGED OR BROKEN. DO NOT ATTEMPT TO RESTERILIZE PEDIATRIC OR NEONATAL TRACHEOSTOMY TUBES. DO NOT EXPOSE TO TEMPERATURES ABOVE 120°F (49°C).

WHEN USING A LASER IN CONJUNCTION WITH THIS OR OTHER TRACHEAL TUBES, CARE MUST BE TAKEN TO AVOID CONTACT OF THE LASER BEAM WITH THE TUBE. SUCH CONTACT, ESPECIALLY IN THE PRESENCE OF OXYGEN-ENRICHED MIXTURES, COULD RESULT IN RAPID COMBUSTION OF THE TUBE WITH HARMFUL THERMAL EFFECTS AND WITH EMISSION OF CORROSIVE AND TOXIC COMBUSTION PRODUCTS, INCLUDING HYDROCHLORIC ACID (HCL).

REPLACEMENT TRACHEOSTOMY TUBES SHOULD BE KEPT AT THE PATIENT'S BEDSIDE.

THE DIFFERENCES IN TECHNIQUE BETWEEN PEDIATRIC AND ADULT TRACHEOTOMIES MUST BE STRICTLY OBSERVED. Reference: 1. Othersen, Jr., H.B., The Pediatric Airway (W.B. Saunders Co.—Philadelphia, 1991) p. 183.

POST FOR USE BY ALL TRACHEOSTOMY CARE PERSONNEL

STERILE. Sterilized with ethylene oxide. For single patient use only.

**CAUTION:** Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

##### **CAUTIONS:**

*The Shiley Tracheostomy Tube and Accessory Products are composed of medical grade plastics which do not contain latex material.*

*Precautions should be taken when discarding any of the tracheostomy tube devices and accessories. Disposal of medical devices should be made in accordance with applicable national regulations for biologically hazardous waste.*

Contents of Packages—Refer to Table 1.

Description—Refer to Tables 2 through 5 For Sizing Dimensions

**NOTE:** The sizing dimensions as noted in Tables 2 through 4 refer to the I.D. as the inside diameter at the smallest location. O.D. is the outside diameter of the cannula (not including cuff features). Length is the distance from the patient side of the neck plate to the patient end of the tube, measured on the centerline.

### Description

The Shiley Pediatric and Neonatal Tracheostomy Tubes have a radiopaque tube constructed using biocompatible polyvinyl chloride with a built-in standard 15mm connector on the proximal end of the tube for direct connection to standard ventilation and anesthesia equipment. Cuffed Pediatric models have a thin wall, high volume, low pressure cuff to minimize tracheal wall damage. When inflated, the cuff conforms to the natural shape of the trachea providing a seal at a low intracuff pressure. The cuff inflation line has a luer valve with an integral pilot balloon to indicate cuff inflation.

### Indications For Use

This device is intended for use in providing tracheal access for airway management.

#### **CAUTION:**

*The Shiley Tracheostomy Tube is classified as a disposable medical device. The manufacturer recommends that the Tracheostomy Tube usage not exceed twenty-nine (29) days. Frequent and routine changes of the Tracheostomy Tube and accessories are recommended and should be evaluated by the attending physician.*

### Instructions

(Additional instructions for cuffed products are given at the end of this section.)

1. The physician should determine whether a Shiley Neonatal or Pediatric tube should be used.

**NOTE:** A Neonatal tube has the same inner and outer diameter as the Pediatric tube of the same size but has a shorter cannula and neck flange. A long Pediatric tube has the same inner and outer diameter as a standard Pediatric tube of the same size, but has a longer cannula.

2. To provide the maximum airway, the largest size tube which will conform to the patient's anatomy should be chosen.
3. Expose the anterior tracheal wall using a standard infant tracheostomy technique. Reference: 1. Othersen, Jr., H.B., *The Pediatric Airway* (W.B. Saunders Co.—Philadelphia, 1991) p. 183.
4. A thin film of water soluble lubricant should be applied to the tube and protruding portion of the obturator to facilitate ease of inserting the tube. Insert the tube with the obturator fully seated in place.
5. After the tube is properly in place, immediately remove the obturator.

6. Attach the neck strap to the tube. Flex the patient's head forward and tie the neck strap. When properly adjusted, you should be able to place one finger between the neck strap and the patient's neck.

#### Additional Instructions For Cuffed Pediatric Tubes

**CAUTIONS:**

*The physician should perform a tracheostomy procedure that avoids damage to the cuff. Avoid sharp edges of cartilage in the incision because they may perforate the cuff during insertion and cause leakage.*

*To ease insertion and to guard against cuff perforation from sharp edges of cartilage, the cuff should be tapered back. This can be accomplished by first inflating the cuff. Then gently move the cuff away from the distal tip of the cannula towards the neck plate as the residual air is removed by deflation. Do not use any sharp instruments such as forceps or hemostats that would damage the cuff when tapering it.*

*This product is composed of soft materials which conform to tracheal tissue for performance and patient comfort. Simple precautions in handling of the Shiley Cuffed Pediatric Tracheostomy Tube during insertion and while in place will facilitate proper function and minimize tears and breaks in the inflation system.*

*Avoid pulling or manipulation of the inflation line, as it is designed to conduct and hold air as part of the cuff inflation system. It is recommended that the inflation line be maintained in a position allowing for patient mobility without placing tension on the line-to-flange junction.*

*Prevent lint or other particulates from entering the luer valve of the pilot balloon. Frequent and routine changes of the Tracheostomy Tube are recommended.*

**WARNINGS:**

**THE CUFFED PEDIATRIC TRACHEOSTOMY TUBE (PDC AND PLC) SHOULD NEVER BE PLUGGED OR CAPPED BECAUSE THERE MAY NOT BE SUFFICIENT AIRWAY FOR THE PATIENT TO BREATHE EVEN WITH THE CUFF COMPLETELY DEFLATED.**

**UNDER NO CIRCUMSTANCES SHOULD MORE THAN 25 mm OF Hg AIR PRESSURE BE USED TO INFLATE THE CUFF. OVER-INFLATION OF THE CUFF MAY CAUSE TRACHEAL DAMAGE AND MAY INHIBIT VENTILATION.**

1. Always test the cuff and inflation system for leakage before inserting the tube. This test can be performed as follows: Inflate the cuff with the volume of air indicated in Table 4. Then either observe for deflation over several minutes or immerse the tube in sterile saline and observe for air leakage. Deflate the cuff prior to insertion.

**NOTE:** Inflation volumes are for test purposes only. Consult the physician or Home Health Care Provider for the appropriate inflation volume when the tube is positioned within the trachea.

2. Inflate the low pressure cuff by injecting air into the luer valve of the inflation line using a syringe. Selection of a cuff inflation and deflation procedure is left to the discretion of the physician. References: 1. Morrison, M.L., *Respiratory Intensive Care Nursing* (Little, Brown & Co.—Boston, 1979) pp. 993-108. 2. Shapiro, B., et al. *Clinical Application for Respiratory Care* (Year Book Medical Publishers—Chicago, 1979) pp. 254-257, 261-280.
3. To deflate the low pressure cuff, withdraw air slowly from luer valve of the inflation line using a syringe.
4. Accumulated secretions above the cuff may need to be suctioned before deflating the cuff with syringe, unless suctioning is contraindicated.

### Tracheostomy Tube Care

**WARNINGS:**

**DO NOT USE SOLUTIONS OTHER THAN THOSE RECOMMENDED BELOW TO CLEAN ANY PART OF THE TRACHEOSTOMY TUBE.**

**DO NOT EXPOSE THE TRACHEOSTOMY TUBE TO ANY CHEMICAL AGENTS OTHER THAN THOSE RECOMMENDED, AS THEY MAY RESULT IN TUBE DAMAGE.**

**DO NOT CLEAN CUFFED PORTION OF CANNULA.**

**DO NOT SOAK CANNULA IN PEROXIDE.**

**CANNULA & OBTURATOR**

1. Clean with hydrogen peroxide (half-strength), household vinegar (half-strength), sterile normal saline, or water and mild detergent.
2. After cleaning, thoroughly rinse with sterile saline to remove any cleaning solution residues.

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Pleasanton, CA U.S.A.
Tyco Healthcare UK LTD