



Comfort. For us, it's personal.

With its new StabilitySelector and ball-and-socket headgear attachments, the **ComfortSelect™** Nasal Mask delivers the kind of secure, personalized fit that your patients demand. Its textured, dual-layer cushion design and adjustable forehead arm help deliver a comfortable fit and guard against mask leaks. It's one more example of how the Comfort Series of Masks helps fit 100% of your patients 100% of the time.

t h e c o m f o r t s e r i e s™

t h e c o m f o r t s e r i e s [™]

ComfortSelect[™]

Adjustable StabilitySelector[™] helps make this Comfort Series mask a comfortable, leak-free fit for your patients – and your business.


The ComfortSelect Nasal Mask is an ideal mask for a wide range of patients. Its adjustable StabilitySelector and soft silicone forehead pad help patients find the right fit to reduce pressure points and mask leaks – without overtightening the headgear. And once you find that perfect fit, the ComfortSelect's ball-and-socket headgear attachments make it easier to keep it – they allow you to remove the mask from the headgear without loosening the headgear straps. Once it's on, even gentle tugging will not result in mask leaks. The fit is that secure. And that makes compliance easier to come by.

We've also made sure that the ComfortSelect is a mask that you can feel secure recommending. It offers many features that have been proven popular in other Respiroics masks, like its dual-layer mask cushion and silent exhalation port. ComfortSelect uses fewer mask sizes to fit the same population as comparable masks, and it has minimal parts to lose or break, both of which are particularly comforting when you have a crowded stockroom to manage.

www.comfortseries.respiroics.com

ordering info

Item	Part #
SMALL MASK WITH HEADGEAR	1007930
SMALL MASK ONLY	1007933
MEDIUM MASK WITH HEADGEAR	1007919
MEDIUM MASK ONLY	1007932
SMALL/WIDE MASK WITH HEADGEAR	1007931
SMALL/WIDE MASK ONLY	1007934
Parts and Accessories	
	Part #
SMALL CUSHION WITH RETAINING RING	1007936
MEDIUM CUSHION WITH RETAINING RING	1007935
SMALL/WIDE CUSHION WITH RETAINING RING	1007937
DELUXE HEADGEAR	1002800
DELUXE HEADGEAR, SMALL	1018755
FOREHEAD SUPPORT AND PAD	1018401
SILICONE FOREHEAD PADS (4/PACK)	1009722
BALL-AND-SOCKET HEADGEAR SWIVEL CLIPS (4/PACK)	1009723
EXHALATION ELBOW KIT	1004395
MASK SWIVEL	7041
OXYGEN ENRICHMENT ATTACHMENT	312710
SIZING GAUGE	1012104



ComfortSelect uses Respiroics' one-size-fits-most **Deluxe Headgear**. This unique headgear with Velcro[®] tabs easily adjusts in eight places for a secure, custom fit. The Deluxe Headgear evenly distributes tension, adds stability and helps to prevent mask leaks. The ComfortSelect and Headgear disassemble easily for cleaning and replacement. The Headgear can be used with most Respiroics masks.


RESPIRONICS[®]
www.respiroics.com

Customer Service: 1-800-345-6443 or 724-387-4000
 Respiroics Europe: +33-1-55-60-19-80
 Respiroics Asia Pacific: +81-3-5800-0724

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician.
 ComfortSelect, Comfort Series, StabilitySelector and Respiroics are trademarks and "People. Products. Programs." is a service mark of Respiroics, Inc. and its affiliates.
 Velcro[®] is a trademark of Velcro Industries. Latex Free • Patent Pending

©2004 Respiroics, Inc. and its affiliates. All rights reserved.

1020940 KW 6/2/04



SLEEP
MANAGEMENT



CHRONIC RESPIRATORY
MANAGEMENT



TOTAL VENTILATION
SOLUTIONS



RESPIRATORY
DRUG DELIVERY



NEONATAL &
INFANT CARE