



We apologize in advance for making your patients oversleep.

With its new durable, lightweight design, the **ComfortFull™ 2** makes it easier than ever for patients who breathe through the mouth to comply with their therapy and get a good night's rest. And that's one more reason you can fully count on the Comfort Series to help you fit 100% of your patients 100% of the time.

t h e c o m f o r t s e r i e s ™

# ComfortFull™ 2

## More Comfort. More Durability.

With the revolutionary ComfortFull™ mask, we improved the lives of patients whose therapy was compromised by mouth breathing. Now, we've perfected the original by creating a more durable full-face mask with increased seal and comfort: ComfortFull 2.

The outer layer of the ComfortFull 2's dual-flap cushion features a soft, polished surface. This improves comfort, and provides a more reliable seal for better therapy delivery. Plus, the redesigned cushion and faceplate are now just one lightweight piece.


Additional product advancements include an entrainment valve with a new split washer for better swiveling of the elbow and a more secure and standardized connection. Plus, the ComfortFull 2 features a StabilitySelector™ that's easier to use.

With a full range of features, ComfortFull 2 provides you with the value of increased patient compliance, less callbacks or complaints, and a wide range of adjustability that allows for fewer sizes and less inventory. Translation: it's the best option available for full-face masks. Which is good for your patients. And great for your business.

**To learn more about the ComfortFull 2 and the entire Comfort Series of Masks, visit [www.comfortseries.respironics.com](http://www.comfortseries.respironics.com).**

ordering info				
Item	small	medium	large	HCPCS
MASK WITH HEADGEAR	1004880	1004872	1004950	A7030 A7035
MASK ONLY	1010868	1010869	1010870	A7030

replacement parts/accessories		
Item	Part Number	HCPCS
COMFORTHEADGEAR	1014333	A7035
DELUXE HEADGEAR	1002800	A7035
DELUXE HEADGEAR, SMALL	1018755	A7035
FOREHEAD SUPPORT AND PAD	1018401	
SILICONE FOREHEAD PADS (4/PACK)	1009722	
ENTRAINMENT VALVE SWIVEL	1004871	
PORT CAP (5/PACK)	1010872	
BALL-AND-SOCKET SWIVEL CLIP (2/PACK)	1009723	
SIZING GAUGE	1012155	



ComfortFull 2 uses Respironics' **ComfortHeadgear™**. This headgear is specifically designed for a sturdier, more stable mask foundation. The ComfortHeadgear can be used with most Respironics masks.



Customer Service: 1-800-345-6443 or 724-387-4000  
 Respironics Europe: +33-1-47-52-30-00  
 Respironics Asia Pacific: +81-3-5280-9611

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician.

ComfortFull, Comfort Series, StabilitySelector, ComfortHeadgear and Respironics are trademarks and "People. Products. Programs." is a registered service mark of Respironics, Inc. and its affiliates.

©2005 Respironics, Inc. and its affiliates. All rights reserved.

BD 12/06/05 MCI 4100406

PN 1030677